

These validated data are a snapshot at a given moment; further updates are always possible.

<u>Species:</u>	<i>Escherichia coli</i>
<u>Group:</u>	K12xB
<u>Strain designation:</u>	HB101
<u>Accession number:</u>	<b>LMBP 684</b>
<u>Deposit date:</u>	01/12/2008
<u>Depositor:</u>	Prof. Dr E. Remaut <sup>1 2</sup> <sup>1</sup> Department for Molecular Biomedical Research, VIB, Ghent, Belgium <sup>2</sup> Department of Biomedical Molecular Biology, Ghent University, Ghent, Belgium ← B. Allet <sup>3</sup> <sup>3</sup> Glaxo Institute for Molecular Biology S.A., Plan-les-ouates, Geneva, Switzerland
<u>Other culture collection numbers:</u>	<a href="#">CGSC 12554</a> <a href="#">ATCC 33694</a> <a href="#">DSM 1607</a>
<u>Containment level:</u>	This strain has been assigned the containment level 'Class 1' following the European Directive 2009/41/EC on the contained use of genetically modified organisms, and its updates (see also the <a href="#">Belgian risk group classification</a> ).
<u>Medium:</u>	LB-Lennox
<u>Selection marker:</u>	streptomycin (25 µg/ml)
<u>Cultivation temperature:</u>	37°C
<u>Original reference:</u>	Boyer et al., J. Mol. Biol. 41 (1969), 459-472 [ <a href="#">PMID: 4896022</a> ]
<u>Related reference:</u>	Sambrook et al. (eds), Molecular Cloning: A Laboratory Manual, Cold Spring Harbor Laboratory Press, NY (1989) [ISSN/ISBN: 0879693096]
<u>Genotype:</u>	<i>F araC14 leuB6(Am) Δ(gpt-proA)62 lacY1 glnX44(AS) galK2(Oc) λ<sup>-</sup> recA13 rpsL20 xylA5 mtl-1 thiE1 [hsdS20]</i> (Source: CGSC 12554)
<u>Phenotype:</u>	Str <sup>R</sup> r <sub>B</sub> <sup>-</sup> m <sub>B</sub> <sup>-</sup>
<u>Properties:</u>	HB101 is a hybrid K12xB strain, used for plasmid propagation. The strain is recA <sup>-</sup> which minimizes recombination events and has a restriction-minus background (hsdS20) that allows better representation when cloning methylated DNA. HB101 is suitable for many molecular biology applications, including cloning with limited amounts of DNA of vectors that do not require α-complementation for blue/white screening.
<u>Restricted use:</u>	<a href="#">BCCM MTA</a>

### **Culture recovery and preservation instructions**

The enclosed culture has been grown overnight to saturation, confirming its viability. BCCM/GeneCorner advises to recover it immediately on receipt before use or storage.

Recovery: subculturing into liquid or solid medium according to the cultivation conditions above

Long-term preservation: lyophilisation  
cryopreservation (at least at -80°C)